Original  
(in German)

Fritsch GmbH Milling & Sizing | Industriestr. 8 | 55743 Idar-Oberstein

20.10.2023

Green Deserts

Herr Dr. Guido Bohmann

████████████████████  
49377 Vechta  
Deutschland

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Green Deserts  
Herr Dr. Guido Bohmann

[REDACTED]  
49377 Vechta  
Deutschland

Your sample: **Desert sand**  
with report **number: M230260**

Dear Dr Bohmann,

Please find enclosed our grinding report on the processing of the material you sent us. Below you will find details of the grinding process we carried out.

Grinding to the desired fineness is easily achievable. We recommend using a PULVERISETTE 5 classic line planetary mill with 4 beaker stations to grind larger sample quantities.

Please do not hesitate to contact us should you have any further queries. Please always quote the above report number.

Kind regards

**FRITSCH GMBH · Grinding and Measuring**

L. Benes  
Laboratory Manager

S. Donauer

This application report is based exclusively on the preparation of the available sample material in the specified quantity.

No legal claims may be derived from this report. Subject to technical changes and errors.

© Fritsch GmbH - Grinding & Measurement - www.fritsch.de - lab@fritsch.de

**Report number M230260 - Result 01****Material sample:** Desert sand – dry grinding**Uüöh/e:****PULVERISETTE 6 classic line planetary single-chamber mill****Equipment:**

**Material Beaker** Agate  
**type** Grinding beaker 500 ml  
**Number of grinding balls** 70 x  
 15 mm Ø grinding balls

**Grinding balls**

Speed	650	revolutions per minute
Feed size	1	mm
Sample quantity	350	g
Duration	10	min
Final fineness	~ 12 % < 2	µm

**Comment:**

It should be possible to grind your sample to  $d_{10} < 2 \mu\text{m}$  without additives, as specified in questionnaire M230260. For our first attempt, a large 500 ml agate grinding jar was filled with 70 x 15 mm B grinding balls. 350 g of desert sand was added and ground at maximum speed for 5 minutes.

We then verified the resulting particle size distribution using static light scattering (ANALYSETTE 22 NeXT Nano). More than 5% of the sample had already been ground to  $< 2 \mu\text{m}$ .

After a total grinding time of 10 minutes, the desired value of 10%  $< 2 \mu\text{m}$  had already been achieved. The sample stuck slightly to the grinding balls and the sides of the beaker. It is likely that the desired fineness of 20%  $< 2 \mu\text{m}$  cannot be achieved (see Result 02).

**Protocol number M230260 - Result 02**
**Material sample:** Desert sand – grinding in suspension

**Titel:**
**PULVERISETTE 6 classic line planetary single-vessel mill**
**Equipment:**

**Material Beaker** Agate  
**type** Grinding beaker 500 ml  
**Number of grinding balls** 70 x  
 15 mm Ø grinding balls

**Grinding balls**

Speed	650	revolutions per minute
Feed size	1	mm
Sample quantity	350	g
Additive	Water	
Amount of additive	115.0	g
Duration	30	min
Final fineness	~ 22% < 2	µm

**Comment:**

In this experiment, your sample was initially dry-ground in the same way as in Result 01. The grinding was carried out in cycles. A 2-minute pause was programmed after every minute of grinding. This was intended to prevent excess pressure building up in the grinding jar, even when grinding in suspension (e.g. with water) (as pressure above approx. 7 bar would force the seal out of position).

After 15 minutes of dry grinding, the fineness of the sample was only marginally better than after 10 minutes of dry grinding in Result 01. We added 115 g of deionised water and continued grinding in suspension for a further 15 minutes.

After a total grinding time of 30 minutes, the fineness was checked. The volume of particles < 2 µm had already risen to 22%.

The duration of wet grinding could certainly be further optimised if, after a short dry grinding period with large balls, a switch were made to 5 mm B grinding balls. This would mean that, with 2000x balls per revolution, the sample would be ground far more effectively than with the 70x 15 mm IB grinding balls.

Direct use of 5 mm IB grinding balls is recommended for a fineness of 300–500 µm.

We sieved the grinding balls (3 mm mesh) and transferred the ground sample suspension into a plastic bottle. Unfortunately, due to time constraints, we are unable to grind 4 kg of sand in our application laboratory.